What is a Bionic

Pancreas?

It's the future of insulin therapy.

No more:

O Carb counting*

Ocrrection factors

Carb ratios

Pre-set basal rates

or any of the other settings that might be overwhelming about other insulin delivery devices.

The iLet only needs one number - your weight.

After startup, the iLet learns about your basal insulin needs, how to correct you to target, and how much insulin you really need for meals — no more guessing!



Based on clinical trial results, adults using the iLet could experience almost 1,000 extra hours of time in range each year.¹



Ready to GO BIONIC?

Scan the code or reach out to your local representatives.

Safety Information

The iLet® Bionic Pancreas System is indicated for use by people with type 1 diabetes 6 years of age and older. The iLet Bionic Pancreas requires prescription by a physician. Refer to the iLet Bionic Pancreas System User Guide at www.betabionics.com/user-guides, or for complete safety information including indications, contraindications, warnings, cautions, compatible devices, compatible drugs and instructions, refer to www. betabionics.com/safety. DO NOT start to use the iLet Bionic Pancreas System without adequate training. Incorrect use may result in over-delivery or under-delivery of insulin, which could lead to hypoglycemia or hyperglycemia.

Medical Disclaimer

This handout is for information only and is not a substitute for medical advice and/or services from a healthcare provider. All personal health care decisions and treatment should be discussed with a healthcare provider who is familiar with your individual needs.

*User must be carb aware.

Clinical Contact

Sales Contact

Beta Bionics

betabionics.com